

Information for all hearing aid users



Users are asked to check the plan of the theatre below, which highlights the seats that are covered by our hearing assistance loop.

Stage

	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
A																	A
B																	B
C																	C
D																	D
E																	E
F																	F
G																	G
	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	
H																	H
J																	J
	1	2	3	4	5	6	7	8		9	10	11	12	13	14	15	16

All hearing aids must be switched to the “T” position to receive the loop signal

Users of digital aids who find they do not get good reception despite doing this and sitting in an area of good coverage are advised to contact their audiologist and check that the hearing aid has been set up for loop systems.

**Our system is monitored and regularly checked for correct operation
If you think there is a problem please let us know**